

PRODUCT DATA SHEET

Published: 2026-07-10

Hose set (Puridest)

for PD 2 R - PD 12 R, PD 2 D - PD 8 D

Part Number: A000138

Features

- Consisting of hoses for water inlet/outlet (1.5 m) and hose clamps/cable ties
- - 1.5 m pressure hose, 13 mm inner diameter, with two hose clamps
- - 1.5 m SIL hose, 10 mm inner diameter, with one cable tie
- - 1.5 m SIL hose, 15 mm inner diameter, with one cable tie



LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser