

PRODUCT DATA SHEET

Published: 2026-07-06

Hose set
NPT 1/2", 10 m
Part Number: A001515

Features

- Includes 2 hoses and 4 sets of fittings



Technical Features

Connection 1 thread (outer)	NPT 1/2
Connection 2 thread (outside)	NPT 1/2
Inner diameters	15 mm
Outer diameter	22 mm
Insulation thickness	9 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser

PRODUCT DATA SHEET

Published: 2026-07-06

Hose set

NPT 1/2", 10 m

Part Number: A001515

CHILLER SIDE

WATER INLET

1 X 1/2" HOSE CARRIER
10m 15-22 FLEXIBLE PIPE
1 X 21-23 CLAMPS
10m ARMAFLEX INSULATION
10m TECNOTHERM TAPE



WATER OUTLET

1 X 1/2" HOSE CARRIER
10m 15-22 FLEXIBLE PIPE
1 X 21-23 CLAMPS
10m ARMAFLEX INSULATION
10m TECNOTHERM TAPE



 Loctite 55

APPLICATION SIDE

WATER INLET

1 X 1/2" HOSE CARRIER
1 X 21-23 CLAMPS

WATER OUTLET

1 X 1/2" HOSE CARRIER
1 X 21-23 CLAMPS

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser