

PRODUCT DATA SHEET

Published: 2026-03-16

Cylindrical insert
suitable for U 20 M and U 2040 M
Part Number: A001722

Features

- 2-chamber construction for constant level in the thermostating chamber
- High temporal and spatial temperature stability thanks to the laminar flow that develops
- Temperature stability with water, at $T = 70\text{ °C}$ and pump stage 5: $<0.01\text{ K}$
- Temperature homogeneity with water, at $T = 70\text{ °C}$ and pump stage 5: $<0.01\text{ K}$
- Two screws, accessible from above, for leveling the cylinder insert
- Please follow the installation instructions!



Technical Features

Outer diameter	184 mm
Inner diameters	175 mm
Immersion depth	300 mm
Material	Stainless steel

Standard accessories

- 1 Bath cover

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser