

PRODUCT DATA SHEET

Published: 2026-10-09

Jet pipe
suitable for U 12 M, U 1245 M
Part Number: A001827

Features

- For improving temperature homogeneity in the bath.
- Excellent temperature stability ($\pm 0.004\text{ }^{\circ}\text{C}$) and uniformity ($\pm 0.003\text{ }^{\circ}\text{C}$). Measured with pump stage 6 and water at $70\text{ }^{\circ}\text{C}$ (according to DIN 12876-2).
- The jet pipe is mounted on the pump outlet and is positioned centrally on the bath floor up to the front of the bath.
- Installation instructions are included.




Technical Features

Material

Stainless steel

Weight

0.20 kg



Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.),
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser