

PRODUCT DATA SHEET

Published: 2026-06-29

Pt 100 Platinum resistance thermometer
 measuring range: -200 °C ... 450 °C, cable: 1.5 m
 Part Number: ETP 065

Features

- Immersion probe, sensor class 1/3 DIN
- With black handle and cable 2 m



Technical Features

Outer diameter	3 mm
Immersion depth	300 mm
Material	CrNi steel (1.4571)
Weight	0.18 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
 Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
 info@lauda.de • www.lauda.de
 WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
 Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
 LAUDA DR. R. WOBSEY Verwaltungs-GmbH
 Sitz Lauda-Königshofen
 Registergericht Mannheim • HRB 560226

Geschäftsführer:
 Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
 Dr. Marc Stricker
 Beirat: Dr. Gerhard Wobser