

PRODUCT DATA SHEET

Published: 2026-09-15

Adapter M30 x 1.5 O to NPT 1/2" O

Part Number: HKA 186-A



Technical Features

Material

Connection 1 thread (outer)

Connection 2 thread (outside)



Weight

Stainless steel

M30 x 1,5

1/2" NPT(F)

0.11 kg



LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser