

PRODUCT DATA SHEET

Published: 2026-08-12

Screw plug for M16 x 1

without union nut

Part Number: HKN 065


Features

- Classic austenitic steel, also commonly referred to as "free-cutting steel".
- Good corrosion resistance, but limited acid resistance and not resistant to intergranular corrosion.



Technical Features

Outer diameter	13.5 mm
Inner diameters	10 mm
Material	CrNiS steel (1.4305)
Length	1 cm
Weight	0.01 kg



LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser