

PRODUCT DATA SHEET

Published: 2026-06-29

Support tray with placement surface
8 steps, with footprint 115 x 270 mm
Part Number: LCZ 0648

Features

- Rising platform made of stainless steel.
- Fixed on one site of the bath.
- Height adjustable in eight steps. Distance between steps: 15 mm
- Not for use with bath cover or gable cover.

Technical Features

Overall dimensions (WxDxH)

297.1 x 118 x 150 mm

Usable floor space width

268.7 mm

Usable floor space depth

115 mm

Material

Stainless steel

Weight

0.60 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser