

PRODUCT DATA SHEET

Published: 2026-06-29

Transparent bath

12 l

Part Number: LCZ 0784

Features

- Transparent bath made from polycarbonate



Technical Features

Overall dimensions (WxDxH)	356 x 331 x 208 mm
Usable floor space width	300 mm
Usable floor space depth	315 mm
Immersion depth	140 mm
Material	Stainless steel
Weight	3.30 kg
Bath opening max. width	300 mm
Bath opening max. depth	315 mm
Bath depth	160 mm
Usable depth	140 mm

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser