

PRODUCT DATA SHEET

Published: 2026-06-29

Automatic refill
with LiBus connection
Part Number: LCZ 9661

Features

- The automatic filling device compensates the loss of evaporated bath liquid, by especially by using baths with a large open bath surface and at higher temperatures.
- When using water as bath liquid, the unit can be connected to a water supply via the automatic filling device. In addition, units with a delivery/suction pump can draw the bath liquid from a storage container (e.g. canister) via the automatic filling device.
- The automatic filling device can be connected to units which are equipped with a LiBus modul and which have a multi-step level measurement.
- The valve is designed for operation down to -10 °C.



Technical Features

Weight

0.68 kg

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser