

PRODUCT DATA SHEET

Published: 2026-09-02


Coolant valve

with LiBus connection, G 3/8" I to G 3/8" O

Part Number: LCZ 9771

Features

- The cooling liquid valve is designed for operation with cooling water only.
- With the cooling liquid valve, the cooling water quantity can be controlled (on/off) by means of LiBus module, by opening the cooling water inlet only when cooling is required.
- For safe operation of the cooling liquid valve with water above 100 °C, the free hose end of the cooling coil must be fixed at the outlet, as a steam blow may occur in the short



Technical Features

Connection 1 thread (inside)

G 3/8"

Connection 2 thread (outside)

G 3/8"

Weight

0.70 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser