

PRODUCT DATA SHEET

Published: 2026-06-26

Sight glass, M38 x 1.5 O to M38 x 1.5 I
-40 °C to 140 °C

Part Number: LWZ 141

Features

- Check the heat transfer liquid by visual inspection, e.g. for degassing status, moisture turbidity, color change



Technical Features

Overall dimensions (WxDxH)

73 x 45 x 182 mm

Operating temperature range

-40 ... 140 °C

Pressure max.

12 bar

Connection 1 thread (inside)

M38 x 1,5

Connection 2 thread (outside)

M38 x 1,5

Weight

0.93 kg

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser