

PRODUCT DATA SHEET

Published: 2026-06-29

Pipe set FC 80 MID / FC 80 C on Integral
Housing like IN 1850 XTW, insulated
Part Number: LWZ 951

Features

- Contains two pipes for hydraulic connection of a FC 80 MID on an Integral such as IN 1850 XTW
- Closed-pore insulating material, wall thickness 25 mm

Technical Features

Operating temperature range

-40 ... 140 °C

Connection 1 thread (inside)


M38 x 1,5

Connection 2 thread (inside)

M38 x 1,5

Material

Stainless steel / Copper



LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser