

## PRODUCT DATA SHEET

Published: 2026-08-26

Pipe set FC 80 MID / FC 80 C on Integral  
Housing like IN 3540 XTW, insulated  
Part Number: LWZ 952

### Features

- Contains two pipes for hydraulic connection of a FC 80 MID or FC 80 C on an Integral such as IN 1850 XTW
- Closed-pore insulating material, wall thickness 25 mm



### Technical Features

Operating temperature range

-40 ... 140 °C

Connection 1 thread (inside)



M38 x 1,5

Connection 2 thread (inside)

M38 x 1,5

Material

Stainless steel / Copper



Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG  
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0  
info@lauda.de • www.lauda.de  
WEEE-Reg.-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen  
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:  
LAUDA DR. R. WOBSEY Verwaltungs-GmbH  
Sitz Lauda-Königshofen  
Registergericht Mannheim • HRB 560226

Geschäftsführer:  
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,  
Dr. Marc Stricker  
Beirat: Dr. Gerhard Wobser