

PRODUCT DATA SHEET

Published: 2026-06-29

Metal hose

M16 x 1, temp. range: 10 ... 400 ° C, 100 cm

Part Number: LZM 041

Features

- Corrugated metal hose for temperature control with different heat transfer liquids
- Conical sealing
- Static bending radius: 50 mm; Dynamic bending radius: 180 mm
- Hose available by the metre



Technical Features

Operating temperature range

10 ... 400 °C

Connection 1 thread (inside)

M16 x 1

Connection 2 thread (inside)

M16 x 1

Length

100 cm

Inner diameters

10 mm

Outer diameter

22.2 mm

Insulation thickness

1.6 mm

Material

Stainless steel

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser