

PRODUCT DATA SHEET

Published: 2026-08-13

EPDM tube

ID 9 mm, insulation 9 mm

Part Number: LZS 019

Features

- EPDM hose with additional insulation for connection to pump nozzles with hose nipples
- Insulation made of highly flexible, closed-cell insulating material based on EPDM.
- The insulating material is weather-, UV- and ozone-resistant and is characterised by high water vapour diffusion resistance and low thermal conductivity.
- EPDM hose not suitable for mineral oils and Ultra 350



Technical Features

Operating temperature range	-35 ... 90 °C
Inner diameters	9 mm
Outer diameter	32 mm
Insulation thickness	9 mm
Material inside	EPDM
Material outside	EPDM

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser