

PRODUCT DATA SHEET

Published: 2026-08-20

Insulating hose

ID 23 mm x WT 10 mm

Part Number: RKJ 009

Features

- Insulating material for cold applications, closed-pore
- The insulating material is weather-, UV- and ozone-resistant and is characterised by high water vapour diffusion resistance and low thermal conductivity.
- Hose available by the metre



Technical Features

Operating temperature range

-50 ... 150 °C

Inner diameters

23 mm

Insulation thickness

10 mm

Material

EPDM

Reserve technical changes

LAUDA DR. R. WOBSEY GMBH & CO. KG
Laudaplatz 1 • 97922 Lauda-Königshofen • DE

T + 49 (0) 9343 503-0
info@lauda.de • www.lauda.de
WEEE-Reg-Nr.: DE 66 42 40 57

Kommanditgesellschaft: Sitz Lauda-Königshofen
Registergericht Mannheim • HRA 560069

Persönlich haftende Gesellschafterin:
LAUDA DR. R. WOBSEY Verwaltungs-GmbH
Sitz Lauda-Königshofen
Registergericht Mannheim • HRB 560226

Geschäftsführer:
Dr. Gunther Wobser (Vors.), Dr. Mario Englert,
Dr. Marc Stricker
Beirat: Dr. Gerhard Wobser